MR Safe

ADDUCTED ARM SCOOP

AM1100

INSTRUCTIONS FOR USE

Issue Date: 03/2026

Revision: C

EC REP

EC Representative

Qualrep Services B.V.
Utrechtseweg 310 Bldg B42
6812 AR Arnhem
Netherlands

Manufacturer

Adept Medical Ltd
2-6 McDonald St,
Morningside
Auckland 1025, New Zealand

adeptmedical.com | 

ADEPT

Medical
Supporting you

Product Registration

Register today to ensure your product is covered by warranty, and for easy access to product services and general support.

www.adeptmedical.com/contact/product-registration

Warranty

Service


Support

Contents

Unboxing	4
Setup	4
Spare Parts	5
Image Table Verification	5
Essential Information	6


ES Spanish	8	NO Norwegian	30
FR French	10	PL Polish	32
IT Italian	12	PT Portuguese	34
DE German	14	RO Romanian	36
NL Dutch	16	RU Russian	38
HR Croatian	18	SR Serbian	40
CS Czech	20	SK Slovak	42
DA Danish	22	SL Slovenian	44
FI Finnish	24	SV Swedish	46
EL Greek	26	TR Turkish	48
HU Hungarian	28		

Unboxing


Setup


1



2

Placement with soft arm pad or similar

Placement with intravenous (IV) or intra arterial (IA) line

Spare Parts

Product Code

Adducted Arm Scoop x 2

AM1100



Component Code

Adducted Arm Support x 1

M6555

Image Table Verification

The product has been risk assessed to operate within the below specifications

Table Angulation: +/- 15° lateral & longitudinal
Table Surface: Flat / Curved
CT/MRI Bore: ≥70cm

Essential Information

Information

Intended Purpose

To support a patient's adducted arm during image-guided medical procedures.

Intended User/Training Requirement

Intended to be used by trained medical professionals.

Patient Target Group

Adult male or female patient requiring diagnostic imaging and image guided therapy.

Contraindications

- Not to be used with patients over 135 kg.
- Not intended for bore size less than 70cm.

Warnings and Cautions

- Ensure IFU is read prior to use.
- Do not lean on installed product.
- Ensure Preparatory Cleaning is conducted.
- Ensure Preparatory Product Check is conducted.
- Ensure soft arm pad is used on support surface.
- Ensure installed product does not catch IV (intravenous) or IA (intra-arterial) line during placement.
- Do not use with patients over 135 kg.

Incident Reporting Guidelines

For product complaints and incidents, please complete form on www.adeptmedical.com/contact/service-and-repairs

Initial Checks & Cleaning

Cleaning

Refer to Disinfection Instructions.


Product Checks

- No visible damage or sharp edges, e.g. cracks.


MRI Safety Checks

Ensure each component below has the 'MR Safe' product sticker and engraving in the below locations:

- Adducted Arm Support:

- Safety Clearance: $\geq 70\text{cm}$

Disinfection

Warnings

- Insufficient cleaning may compromise disinfection process and lead to transmission of pathogens.
- Do not use any disinfectant products or abrasive/corrosive agents which are not on the approved cleaners list.
- Always read manufacturer's instructions and consult the manufacturer's MSDS for cleaning and disinfectant products.
- Do not use any cleaning procedure which exceeds $65^{\circ}\text{C}/149^{\circ}\text{F}$.

Limitations on Processing

Discontinue use if:

- Any cracks, breakages or tears are present.

Disinfection Instructions

Initial Treatment at the Point of Use

It is important to clean the product once it is removed from the packaging and after each procedure. Ensure all areas including slots are thoroughly cleaned after use to remove all contaminant build-up that may be present following a procedure.

Preparation Before Cleaning

- Remove Adducted Arm Scoop from beneath the mattress
- www.adeptmedical.com/imaging-table-accessories/adducted-arm-scoop/training?group=cleaning-instructions

Cleaning: Manual

- If any visible residue is present, rinse under water.
- A soft bristle brush can be used.
- Ensure the device is completely dried prior to disinfection.

Disinfection

Refer to the Approved Cleaners List for selection of an appropriate disinfecting agent.

- With a cloth dampened by an approved disinfectant, wipe all surfaces being careful to wipe within slots.

Approved Cleaners List

Approved Disinfectant by Brands

Liquids:

ORION Laboratories - 70% Isopropyl Alcohol
Jaychem Industries - 2% Chlorhexidine / 70% Alcohol
Betadine - Povidone Iodine (7.5%)
BODE Chemie GmbH - Kohrsolin FF
STERIS - Coverage Spray TB
CaviCide - Metrex Research

Dissolvable:

Du Pont - Rely+On Virkon Tablets

Surface Wipes:

PDI Sani-Cloth Bleach
PDI Sani-Cloth AF3
Reynard Health Supplies - Surface Disinfectant Wipes
Metrex Research - CaviWipes
Clinell Wipes - Universal (Green)
Clorox Hydrogen Peroxide Wipes

Approved Active Disinfectant Ingredients

Acid Based:

$\leq 10\%$ Malic Acid CAS 6915-15-7
 $\leq 6\%$ Sulfamic acid CAS 5329-14-6

Alcohol Based:

$\leq 5\%$ 2-Butoxyethanol CAS 111-76-2
 $\leq 10\%$ Butyldiglycol CAS 112-34-5
 $\leq 70\%$ Isopropyl Alcohol (Propan-2-ol) CAS 67-63-0
 $\leq 10\%$ Tridecanol CAS 69011-36-5
 $< 10\%$ Alcohols, C12-14, ethoxylated CAS 68439-50-9
 $\leq 70\%$ Denatured Ethanol CAS 64-17-5

Sulphate Based:

$\leq 0.1\%$ PHMB CAS 27083-27-8
 $\leq 55\%$ Potassium Peroxymonosulfate CAS 70693-62-8
 $\leq 3\%$ Potassium Persulfate CAS 7727-21-1

Ammonium and Chloride Based:

$\leq 5\%$ Benzalkonium Chloride CAS 68424-85-1
 $\leq 0.28\%$ Benzethonium Chloride CAS 121-54-0
 $\leq 10\%$ Benzyl-C23-18-Alkyl-dimethyl Ammonium Chloride CAS 8001-54-5
 $\leq 2\%$ Chlorhexidine CAS 55-56-1
 $\leq 10\%$ Didecyl Dimethyl Ammonium chloride CAS 7173-51-5
 $\leq 0.5\%$ Quaternary Ammonium Compounds CAS 68956-79-6
 $< 5\%$ Benzyl-C 12-18 alkyldimethylammonium chlorides CAS 63891-01-5

Others:

$\leq 10\%$ Glutaral CAS 111-30-8
 $\leq 0.63\%$ Sodium Hypochlorite CAS 7681-52-9
 $\leq 7.5\%$ Povidone Iodine CAS 25655-41-8
0.5% Disodium Cocoampho Dipropionate CAS68604-71-7
 $\leq 1.4\%$ Hydrogen Peroxide CAS 7722-84-1

Access the latest Approved Cleaners list via www.adeptmedical.com/resources

Drying

All components should be dried thoroughly before use.

Maintenance, Inspection and Testing

All components:

- Ensure no visible damage, including sharp edges or cracks.

Packaging

It is not required to package the device following disinfection.

Sterilisation

This device should not be subjected to sterilisation processes.

Storage

Once disinfection is complete and all components are dry the device should be stored in a dry environment in which the temperature will not exceed $65^{\circ}\text{C}/149^{\circ}\text{F}$.

Disposal

Used product is a biohazard, decontaminate according to instruction provided in the Disinfection Instructions and reach out to the manufacturer on adeptmedical@adept.co.nz for material specification, if required. Dispose according to internal clinic policy taking into consideration local regulations.

Serious Incidents

Any serious incident which occurs in relation to the device should be reported immediately to:

- The Competent Authority of the applicable Member State.
- The Manufacturer using the Serious Incident Reporting Form found on the Adept Medical Website: www.adeptmedical.com/contact/serious-incident-reporting-form

The instructions provided above have been validated by the manufacturer of the medical device as being capable of preparing a medical device for reuse. It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process.

Visit the Adept Medical Website for detailed visual instructions related to this device. www.adeptmedical.com

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