CLINICIAN ARM SUPPORT

AM0900

INSTRUCTIONS FOR USE

Issue Date: 03/2026

Revision: C

EC REP

EC Representative

Qualrep Services B.V.
Utrechtseweg 310 Bldg B42
6812 AR Arnhem
Netherlands

Manufacturer

Adept Medical Ltd
2-6 McDonald St,
Morningside
Auckland 1025,
New Zealand

adeptmedical.com | 

ADEPT

Medical
Supporting you

Product Registration

Register today to ensure your product is covered by warranty, and for easy access to product services and general support.

www.adeptmedical.com/contact/product-registration

Warranty

Service


Support

Contents

Unboxing	4
Assembly	4
Setup	5
Spare Parts	7
Image Table Verification	7
Essential Information	8


ES Spanish	10
FR French	12
IT Italian	14
DE German	16
NL Dutch	18

Unboxing

	(A)	Arm Support Pad
	(B)	Forearm Assembly
	(C)	Tube Assembly
	(D)	Base Board Assembly

Assembly


5


Setup


Table Setup


Height Adjustment


Angle Adjustment (+/-20°)


Armrest Placement


Stability Range Over the Table

Sterility

Friction Adjustment

Spare Parts


Product Code

	Clinician Arm Support	AM0900
	Arm Support Pad	AM0171

Component Codes

	D	Base Board Assembly	M9100
---	----------	---------------------	-------

Tube Assembly

Forearm Assembly (M9120)



Image Table Verification

The product has been risk assessed to operate within the below specifications

Table Angulation: +/- 20° longitudinal & 10° lateral
Table Surface: Flat

Essential Information

Information

Intended Purpose

To support the Clinicians arm during endoscope guided surgery.

Intended User/Training Requirement

Intended to be used by trained medical professionals.

Patient Target Group

Patients undergoing endoscopic guided surgery.

Contraindications

- Not to be used by clinicians weighing over 135 kg.
- Not to be used with patients weighing under 45 kg.

Warnings and Cautions

- Ensure IFU is read prior to use.
- Ensure Preparatory Cleaning is conducted.
- Ensure Preparatory Product Check is conducted.
- Restless patients could dislodge Clinician Arm Support and need to be monitored.

Incident Reporting Guidelines

For product complaints and incidents, please complete form on: www.adeptmedical.com/contact/service-and-repairs

Initial Checks & Cleaning

Cleaning

Refer to Disinfection Instructions.

Product Checks

- No visible damage or sharp edges.
- Ensure the Tube Assembly is free to move up, down and rotate when the lever is in the unlocked position.
- Ensure the Tube Assembly cannot move in any direction when the lever is in the locked position.
- Ensure all rotating joints can be locked/ unlocked using thumbscrew.
- Arm Support Pad: Ensure no visible damage or tearing.

Disinfection

Warnings

- Insufficient cleaning may compromise disinfection process and lead to transmission of pathogens.
- Do not use any disinfectant products or abrasive/corrosive agents which are not on the approved cleaners list.
- Always read manufacturer's instructions and consult the manufacturers MSDS for cleaning and disinfectant products.
- Some disinfectants may cause slight discolouration to the soft blue material used on some components within the product range. This will not affect the strength and the product will remain fit for purpose.

Limitations on Processing

Discontinue use if:

- Any cracks, tears, or breakages are present.
- Any sharp edges develop.
- If the locking mechanism no longer holds the Tube Assembly in place when locked.

Disinfection Instructions

Initial treatment at the point of use

It is important to clean the product once it is removed from the packaging and after each procedure. Ensure all areas including joints, clips, sockets, brackets and levers are thoroughly cleaned after use to remove all contaminant build-up that may be present following a procedure.

Preparation before cleaning

- Remove the Arm Support Pad
- Remove the Base Board from beneath the mattress

www.adeptmedical.com/imaging-table-accessories/clinician-arm-support/training?group=cleaning-instructions

Cleaning: Manual

- If any visible residue is present, rinse under water, but avoid submerging the device.
- A soft bristle brush can be used.
- Ensure the device is completely dried prior to disinfection.

Disinfection

Refer to the Approved Cleaners List for selection of an appropriate disinfecting agent.

- With a cloth dampened by an approved disinfectant, wipe all surfaces.
- With the lever in the unlocked position, move the Tube Assembly up and down while wiping to ensure the entire surface has been disinfected.
- Ensure the Forearm Assembly is fully expanded to ensure all surfaces are disinfected.

Approved Cleaners List

Approved Disinfectant by Brands

Liquids:

ORION Laboratories - 70% Isopropyl Alcohol
Jaychem Industries - 2% Chlorhexidine / 70% Alcohol
Betadine - Povidone Iodine (7.5%)
BODE Chemie GmbH - Kohrsolin FF
STERIS - Coverage Spray TB
CaviCide - Metrex Research

Dissolvable:

Du Pont - Rely+On Virkon Tablets

Surface Wipes:

PDI Sani-Cloth Bleach
PDI Sani-Cloth AF3
Reynard Health Supplies - Surface Disinfectant Wipes
Metrex Research - CaviWipes
Clinell Wipes - Universal (Green)
Clorox Hydrogen Peroxide Wipes

Approved Active Disinfectant Ingredients

Acid Based:

≤10% Malic Acid CAS 6915-15-7
≤6% Sulfamic acid CAS 5329-14-6

Alcohol Based:

≤5% 2-Butoxyethanol CAS 111-76-2
≤10% Butyldiglycol CAS 112-34-5
≤70% Isopropyl Alcohol (Propan-2-ol) CAS 67-63-0
≤10% Tridecanol CAS 69011-36-5
≤10% Alcohols, C12-14, ethoxylated CAS 68439-50-9
≤70% Denatured Ethanol CAS 64-17-5

Sulphate Based:

≤0.1% PHMB CAS 27083-27-8
≤55% Potassium Peroxymonosulfate CAS 70693-62-8
≤3% Potassium Persulfate CAS 7727-21-1

Ammonium and Chloride Based:

≤5% Benzalkonium Chloride CAS 68424-85-1
≤0.28% Benzethonium Chloride CAS 121-54-0
≤10% Benzyl-C23-18-Alkyl-dimethyl Ammonium Chloride CAS 8001-54-5
≤2% Chlorhexidine CAS 55-56-1
≤10% Didecyl Dimethyl Ammonium chloride CAS 7173-51-5
≤0.5% Quaternary Ammonium Compounds CAS 68956-79-6
≤5% Benzyl-C 12-18 alkyl dimethyl ammonium chlorides CAS 63891-01-5

Others:

≤10% Glutaral CAS 111-30-8
≤0.63% Sodium Hypochlorite CAS 7681-52-9
≤7.5% Povidone Iodine CAS 25655-41-8
0.5% Disodium Cocoampho Dipropionate CAS68604-71-7
≤1.4% Hydrogen Peroxide CAS 7722-84-1

Access the latest Approved Cleaners list via www.adeptmedical.com/resources

Drying

All components should be dried thoroughly before use.

Maintenance, Inspection and Testing

- Ensure there are no visible damage or sharp edges.
- Ensure the Tube Assembly is free to move up and down when the lever is in the unlocked position.
- Ensure the Tube Assembly cannot move when the lever is in the locked position.
- Arm Support Pad: Ensure no visible damage or tearing.

Packaging

It is not required to package the device following disinfection.

Sterilisation

This device should not be subjected to sterilisation processes.

Storage

Once disinfection is complete and all components are dry the device should be stored in a dry environment.

Disposal

Used product is a biohazard, decontaminate according to instruction provided in the Disinfection Instructions and reach out to the manufacturer on adeptmedical@adept.co.nz for material specification, if required. Dispose according to internal clinic policy taking into consideration local regulations.

Serious Incidents

Any serious incident which occurs in relation to the device should be reported immediately to:

- The Competent Authority of the applicable Member State.
- The Manufacturer using the Serious Incident Reporting Form found on the Adept Medical Website: www.adeptmedical.com/contact/serious-incident-reporting-form

The instructions provided above have been validated by the manufacturer of the medical device as being capable of preparing a medical device for reuse. It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process.

Visit the Adept Medical Website for detailed visual instructions related to this device. www.adeptmedical.com

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