HEAD IMMOBILISER

AM0800 • AM0820

INSTRUCTIONS FOR USE

Issue Date: 03/2026

Revision: D

EC Representative
Qualrep Services B.V.
Utrechtseweg 310 Bldg B42
6812 AR Arnhem
Netherlands

Manufacturer
Adept Medical Ltd
2-6 McDonald St,
Morningside
Auckland 1025,
New Zealand

adeptmedical.com |

ADEPT

Medical
Supporting you

Product Registration

Register today to ensure your product is covered by warranty, and for easy access to product services and general support.

www.adeptmedical.com/contact/product-registration

Warranty

Service


Support


Contents

- Unboxing 4
- Assembly 4
- Product Use 5
- Setup 7
- Product Configurations 11
- Spare Parts 12
- Image Table Verification 13
- Essential Information 14

ES Spanish 16	NO Norwegian 38
FR French 18	PL Polish 40
IT Italian 20	PT Portuguese 42
DE German 22	RO Romanian 44
NL Dutch 24	RU Russian 46
HR Croatian 26	SR Serbian 48
CS Czech 28	SK Slovak 50
DA Danish 30	SL Slovenian 52
FI Finnish 32	SV Swedish 54
EL Greek 34	TR Turkish 56
HU Hungarian 36	

Unboxing

	(A)	Hood Assembly x 1
	(B)	Base Board Assembly x 1
 x10	(C)	Forehead Straps x 10


Head Immobiliser (AM0800)

(A) (B) (C)


Hood Assembly + Forehead Straps (AM0820)

(A) (C)

Assembly


Hood Assembly Removal


Product Use

Device Position for Different Tables


Placement on a Curved Table

Use foam cushion, or similar

Device Adjustments

Height Adjustment

Rotational Adjustment

Neck Flexion/Extension

Neck Flexion


Neck Extension (Optional Pillow)

Draping (Optional)


Valve Operation

Connect Hood Port to Central Hospital Vacuum via Yankauer Suction Tube (or similar)

Disconnecting Hood


Setup


1


2 Manage hair


3


4



Beads (filling)

5

6



7


Adjust the device


8

Evacuate hood


Patient Release


Product Configurations

STARSystem (AM0160) + Adducted Arm Scoop (AM1100)
Not Included

Retrograde IR Platform (AM0400) + Adducted Arm Scoop (AM1100)
Not Included

Spare Parts

Product Codes


	Head Immobiliser	AM0800
	Hood Assembly + Forehead Straps x10	AM0820

Component Codes


	Base Board Assembly x1	M8000
	Base Board	M8100
	Hood Support	M8200
	Adjustable Assembly	M8300
	Forehead Straps x 10	AM0830

Hood Assembly Beads

If any beads appear outside the hood, DO NOT USE.



Disassembly

Maintenance instructions:

www.adeptmedical.com/imaging-table-accessories/head-immobiliser/training?group=maintenance-instructions

Image Table Verification

The product has been risk assessed to operate within the below specifications

Table Angulation: +/- 15° lateral & longitudinal
Table Surface: Flat / Curved
CT Bore: ≥70cm

Essential Information

Information

Intended Purpose

To support, stabilise and position a patient's head during image-guided (CT and Angiography) medical procedures.

Intended User/Training Requirement

Intended to be used by trained medical professionals.

Patient Target Group

The intended patient is an Adult male or female requiring diagnostic imaging and image guided therapy.

Contraindications

- Not to be used with patients over 135 kg.
- Not to be used for cases longer than 6 hours duration.
- Not for paediatric use.
- Not to be used on patients with Parietal thinning.
- Not to be used on patient with focal calvaria thinning and other skull abnormality.
- Not to be used on patient with severe cervical spine disease (osteoarthritis or rheumatoid arthritis or spondylarthritis).
- Not to be used on patients that have had craniectomy or craniotomy performed in the last 2 years.
- Not to be used on patient with recent surgery to the scalp, may affect wound healing.
- Not to be used on patient with recent traumatic injury to the head.

Warnings and Cautions

- Ensure IFU is read prior to use.
- Ensure Preparatory Cleaning is conducted.
- Ensure Preparatory Product Check is conducted.
- Do not use with patients over 135 kg.
- Vacuum only, operating vacuum:

Minimum
-50 kPa
-375 mmHG
-15 inHG

- Check patient for pressure related injury after 2 hours.
- 2 person set up.
- Ensure head mass is supported if unlocking cam lever and adjusting position.
- Not for paediatric use.
- Remove hair clips, jewellery and eyewear.
- Ensure long hair is flat against the back.

Incident Reporting Guidelines

For product complaints and incidents, please complete form on www.adeptmedical.com/contact/service-and-repairs

Initial Checks & Cleaning

Cleaning

Refer to Disinfection Instructions.

Product Checks

- No break in Carbon Fibre surface.
- No visible damage or sharp edges, e.g. cracks.
- Lever locked: Ensure Hood Support does not move up/down or left/right when pressure is applied.
- Lever unlocked: Ensure Hood Support can be moved through full range of positioning options.
- Ensure that Hood is locked on to Hood Support.
- Ensure that the Hood has no tears in surface, no beads are visible.
- Ensure Hood sets rigid after 10 seconds of suction.

Disinfection

Warnings

- Insufficient cleaning may compromise disinfection process and lead to transmission of pathogens.
- Do not use any disinfectant products or abrasive/corrosive agents which are not on the approved cleaners list.
- Always read manufacturer's instructions and consult the manufacturer's MSDS for cleaning and disinfectant products.
- Do not use any cleaning procedure which exceeds 65°C/149°F.
- Some disinfectants may cause slight discolouration to the soft blue material used on some components within the product range. This will not affect the strength and the product will remain fit for purpose.

Limitations on Processing

Discontinue use if:

- Any cracks, breakages or tears are present.
- Any exposed carbon fibre is visible.
- Any beads are visible through tear in hood

Disinfection Instructions

Initial Treatment at the Point of Use

It is important to clean the product once it is removed from the packaging and after each procedure. Ensure all areas including joints, sockets and levers are thoroughly cleaned after use to remove all contaminant build-up that may be present following a procedure.

Preparation Before Cleaning

- Remove Head Immobiliser from beneath the mattress
 - Separate Hood Assembly from Base Board Assembly
- www.adeptmedical.com/imaging-table-accessories/head-immobiliser/training?group=cleaning-instructions

Cleaning: Manual

- If any visible residue is present, rinse under water, avoid submerging the hood assembly and keep valve dry.
- A soft bristle brush can be used.
- Ensure the device is completely dried prior to disinfection.

Disinfection

Refer to the Approved Cleaners List for selection of an appropriate disinfecting agent.

Base Assembly

- With a cloth dampened by an approved disinfectant, wipe all surfaces of Base Board, Hood Support and Adjustable Assembly.
- Adjust height of Hood Support to access all Base Board surface.

Hood Assembly

- With a cloth dampened by an approved disinfectant, wipe all surfaces of Hood Assembly.

Approved Cleaners List

Approved Disinfectant by Brands

Liquids:

ORION Laboratories - 70% Isopropyl Alcohol
 Jaychem Industries - 2% Chlorhexidine / 70% Alcohol
 Betadine - Povidone Iodine (7.5%)
 BODE Chemie GmbH - Kohrsolin FF
 STERIS - Coverage Spray TB
 CaviCide - Metrex Research

Dissolvable:

Du Pont - Rely+On Virkon Tablets

Surface Wipes:

PDI Sani-Cloth Bleach
 PDI Sani-Cloth AF3
 Reynard Health Supplies - Surface Disinfectant Wipes
 Metrex Research - CaviWipes
 Clinell Wipes - Universal (Green)
 Clorox Hydrogen Peroxide Wipes

Approved Active Disinfectant Ingredients

Acid Based:

≤10% Malic Acid CAS 6915-15-7
 ≤6% Sulfamic acid CAS 5329-14-6

Alcohol Based:

≤5% 2-Butoxyethanol CAS 111-76-2
 ≤10% Butyldiglycol CAS 112-34-5
 ≤70% Isopropyl Alcohol (Propan-2-ol) CAS 67-63-0
 ≤10% Tridecanol CAS 69011-36-5
 <10% Alcohols, C12-14, ethoxylated CAS 68439-50-9
 ≤70% Denatured Ethanol CAS 64-17-5

Sulphate Based:

≤0.1% PHMB CAS 27083-27-8
 ≤55% Potassium Peroxymonosulfate CAS 70693-62-8
 ≤3% Potassium Persulfate CAS 7727-21-1

Ammonium and Chloride Based:

≤5% Benzalkonium Chloride CAS 68424-85-1
 ≤0.28% Benzethonium Chloride CAS 121-54-0
 ≤10% Benzyl-C23-18-Alkyl-dimethyl Ammonium Chloride CAS 8001-54-5
 ≤2% Chlorhexidine CAS 55-56-1
 ≤10% Didecyl Dimethyl Ammonium chloride CAS 7173-51-5
 ≤0.5% Quaternary Ammonium Compounds CAS 68956-79-6
 <5% Benzyl-C 12-18 alkyl dimethyl ammonium chlorides CAS 63891-01-5

Others:

≤10% Glutaral CAS 111-30-8
 ≤0.63% Sodium Hypochlorite CAS 7681-52-9
 ≤7.5% Povidone Iodine CAS 25655-41-8
 0.5% Disodium Cocoampho Dipropionate CAS68604-71-7
 ≤1.4% Hydrogen Peroxide CAS 7722-84-1

Access the latest Approved Cleaners list via www.adeptmedical.com/resources

Drying

All components should be dried thoroughly before use.

Maintenance, Inspection and Testing

All components:

- Ensure no exposed Carbon Fibre is visible.
- Ensure no visible damage or sharp edges, e.g. cracks.

For the Base assembly:

- Unlock lever and ensure all joints are free to adjust / rotate.

For the Hood assembly:

- Ensure there are no tears in the hood or visible beads.

Packaging

It is not required to package the device following disinfection.

Sterilisation

This device should not be subjected to sterilisation processes.

Storage

Once disinfection is complete and all components are dry the device should be stored in a dry environment in which the temperature will not exceed 65°C/149°F.

Disposal

Used product is a biohazard, decontaminate according to instruction provided in the Disinfection Instructions and reach out to the manufacturer on adeptmedical@adept.co.nz for material specification, if required. Dispose according to internal clinic policy taking into consideration local regulations.

Serious Incidents

Any serious incident which occurs in relation to the device should be reported immediately to:

- The Competent Authority of the applicable Member State.
- The Manufacturer using the Serious Incident Reporting Form found on the Adept Medical Website: www.adeptmedical.com/contact/serious-incident-reporting-form

The instructions provided above have been validated by the manufacturer of the medical device as being capable of preparing a medical device for reuse. It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process.

Visit the Adept Medical Website for detailed visual instructions related to this device. www.adeptmedical.com

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