IR SYSTEM

AM0180 • AM0500

INSTRUCTIONS FOR USE

Issue Date: 03/2026

Revision: H

EC Representative
Qualrep Services B.V.
Utrechtseweg 310 Bldg B42
6812 AR Arnhem
Netherlands

Manufacturer
Adept Medical Ltd
2-6 McDonald St,
Morningside
Auckland 1025, New Zealand

adeptmedical.com |

ADEPT

Medical
Supporting you

Product Registration

Register today to ensure your product is covered by warranty, and for easy access to product services and general support.

www.adeptmedical.com/contact/product-registration

Warranty

Service



Support

Contents

Unboxing	4
STARBoard	4
Assembly	4
Setup	5
IR Shield	6
Setup	6
Spare Parts	7
Image Table Verification	8
Essential Information	10

ES Spanish	12	NO Norwegian	34
FR French	14	PL Polish	36
IT Italian	16	PT Portuguese	38
DE German	18	RO Romanian	40
NL Dutch	20	RU Russian	42
HR Croatian	22	SR Serbian	44
CS Czech	24	SK Slovak	46
DA Danish	26	SL Slovenian	48
FI Finnish	28	SV Swedish	50
EL Greek	30	TR Turkish	52
HU Hungarian	32		

Unboxing

	(A)	STARBoard x1
	(B)	Extension Tray x1
	(C)	Arm Support Pad x1
	(D)	X-Ray Shield x1
	(E)	Daggerboard x1

AM0180 (STARBoard + Extension Tray)

(A) (B) (C)


AM0500 (IR Shield)

(D) (E)



STARBoard



Assembly


Check friction joints and adjust if necessary


Setup


Draping (Optional)

Taping (Optional)


Hand Positioning Options


IR Shield

Setup


1

2

3


Spare Parts

Product Codes

	STARBoard + Extension Tray	AM0180
	IR Shield	AM0500
	Arm Support Pad	AM0171


Component Codes

	Extension Tray	M3334
	X-Ray Shield	AM0630
	Daggerboard	M3677


STARBoard


1


2


3


4

5

6

Maintenance instructions:
www.adeptmedical.com/imaging-table-accessories/ir-system/training?group=maintenance-instructions

STARBoard Component Codes



Image Table Verification

The product has been risk assessed to operate within the below specifications

Table Angulation: +/- 15° lateral & longitudinal
Table Surface: Flat

Essential Information

Information

Intended Purpose

- STARBoard: To support a supine patient's arm, facilitating hyperextension for radial access.
- Extension Tray: To provide a transitional work surface from the radial access site to trolley.
- IR Shield: To provide scatter radiation protection for clinician.

Intended User/Training Requirement

Intended to be used by trained medical professionals.

Patient Target Group

Adult male or female patient requiring diagnostic imaging and image guided therapy.

Contraindications

- Not to be used with patients over 135 kg.
- Not for paediatric use.

Warnings and Cautions

- Ensure IFU is read prior to use.
- Ensure preparatory Cleaning is conducted.
- Ensure Preparatory Product Check is conducted.
- Do not use with patients over 135 kg.
- To provide restraint for restless patient (refer to page 6).
- Not for paediatric use.

Incident Reporting Guidelines

For product complaints and incidents, please complete form on: www.adeptmedical.com/contact/service-and-repairs

Initial Checks & Cleaning

Cleaning

Refer to Disinfection Instructions.

Product Checks

- No exposed Carbon Fibre.
- No break in Carbon Fibre surface.
- No visible damage or sharp edges, e.g. cracks.

STARBoard

- Ensure all pivoting joints can rotate with the appropriate amount of friction and adjust if necessary (refer to page 4).
- Correct placement (refer to page 5).
- Patient can access the handle to grip (refer to page 5).

IR Shield

- Ensure correct daggerboard placement (refer to page 6)
- Ensure the X-Ray Shield is correctly placed in the daggerboard (refer to page 6).

Arm Support Pad

- Ensure no damage or tearing.

Disinfection

Warnings

- Insufficient cleaning may compromise disinfection process and lead to transmission of pathogens.
- Do not use any disinfectant products or abrasive/corrosive agents which are not on the approved cleaners list.
- Always read manufacturer's instructions and consult the manufacturers MSDS for cleaning and disinfectant products.
- Do not use any cleaning procedure which exceeds 65 °C/149 °F.
- Some disinfectants may cause slight discolouration to the soft blue material used on some components within the product range. This will not affect the strength and the product will remain fit for purpose.

Limitations on Processing

Discontinue use if:

- Any cracks, breakages, or tears are present.
- Any exposed carbon fibre is visible.

Disinfection Instructions

Initial treatment at the point of use

It is important to clean the product once it is removed from the packaging and after each procedure. Ensure all areas including joints, clips, sockets, brackets and levers are thoroughly cleaned after use to remove all contaminant build-up that may be present following a procedure.

Preparation before cleaning

- Extension Tray: Detach the Extension Tray from the STARBoard.
- STARBoard: Remove STARBoard from beneath the mattress.
- Arm Support Pad: Remove Arm Support Pad from STARBoard.
- IR Shield: Remove from Daggerboard.
- Daggerboard: Remove from beneath the mattress.

www.adeptmedical.com/imaging-table-accessories/ir-system/training?group=cleaning-instructions

Cleaning: Manual

- If any visible residue is present, rinse under water, but avoid submerging the device.
- A soft bristle brush can be used.
- Ensure the device is completely dried prior to disinfection.

Disinfection

Refer to the Approved Cleaners List for selection of an appropriate disinfecting agent.

STARBoard

- With a cloth dampened by an approved disinfectant, wipe all surfaces, Armrest, Handle, Wrist Pivot, and adjust device to wipe entire Slider Arm.

Extension Tray

- Wipe all surfaces and the bracket area with a cloth dampened by an approved disinfectant.

IR Shield

- With a cloth dampened by an approved disinfectant, wipe all surfaces including inside the slot on the Daggerboard.

Arm Support Pad

- With a cloth dampened by an approved disinfectant, wipe all surfaces.

Approved Cleaners List

Approved Disinfectant by Brands

Liquids:

ORION Laboratories - 70% Isopropyl Alcohol
Jaychem Industries - 2% Chlorhexidine / 70% Alcohol
Betadine - Povidone Iodine (7.5%)
BODE Chemie GmbH - Kohrsolin FF
STERIS - Coverage Spray TB
CaviCide - Metrex Research

Dissolvable:

Du Pont - Rely+On Virkon Tablets

Surface Wipes:

PDI Sani-Cloth Bleach
PDI Sani-Cloth AF3
Reynard Health Supplies - Surface Disinfectant Wipes
Metrex Research - CaviWipes
Clinell Wipes - Universal (Green)
Clorox Hydrogen Peroxide Wipes

Approved Active Disinfectant Ingredients

Acid Based:

≤10% Malic Acid CAS 6915-15-7
≤6% Sulfamic acid CAS 5329-14-6

Alcohol Based:

≤5% 2-Butoxyethanol CAS 111-76-2
≤10% Butyldiglycol CAS 112-34-5
≤70% Isopropyl Alcohol (Propan-2-ol) CAS 67-63-0
≤10% Tridecanol CAS 69011-36-5
<10% Alcohols, C12-14, ethoxylated CAS 68439-50-9
≤70% Denatured Ethanol CAS 64-17-5

Sulphate Based:

≤0.1% PHMB CAS 27083-27-8
≤55% Potassium Peroxymonosulfate CAS 70693-62-8
≤3% Potassium Persulfate CAS 7727-21-1

Ammonium and Chloride Based:

≤5% Benzalkonium Chloride CAS 68424-85-1
≤0.28% Benzethonium Chloride CAS 121-54-0
≤10% Benzyl-C23-18-Alkyl-dimethyl Ammonium Chloride CAS 8001-54-5
≤2% Chlorhexidine CAS 55-56-1
≤10% Didecyl Dimethyl Ammonium chloride CAS 7173-51-5
≤0.5% Quaternary Ammonium Compounds CAS 68956-79-6
<5% Benzyl-C 12-18 alkyl dimethyl ammonium chlorides CAS 63891-01-5

Others:

≤10% Glutaral CAS 111-30-8
≤0.63% Sodium Hypochlorite CAS 7681-52-9
≤7.5% Povidone Iodine CAS 25655-41-8
0.5% Disodium Cocoampho Dipropionate CAS68604-71-7
≤1.4% Hydrogen Peroxide CAS 7722-84-1

Access the latest Approved Cleaners list via www.adeptmedical.com/resources

Drying

All components should be dried thoroughly before use.

Maintenance, Inspection and Testing

All components:

- Ensure no visible damage or sharp edges, e.g. cracks.

STARBoard, X-Ray Shield:

- Ensure no exposed Carbon Fibre is visible.

STARBoard:

- Ensure all pivoting joints can rotate with the appropriate amount of friction and adjust if necessary.

Extension Tray:

- Ensure the Extension Tray is correctly attached and secure on the STARBoard.

Packaging

It is not required to package the device following disinfection.

Sterilisation

This device should not be subjected to sterilisation processes.

Storage

Once disinfection is complete and all components are dry the device should be stored in a dry environment.

Disposal

Used product is a biohazard, decontaminate according to instruction provided in the Disinfection Instructions and reach out to the manufacturer on adeptmedical@adept.co.nz for material specification, if required. Dispose according to internal clinic policy taking into consideration local regulations.

Serious Incidents

Any serious incident which occurs in relation to the device should be reported immediately to:

- The Competent Authority of the applicable Member State.
- The Manufacturer using the Serious Incident Reporting Form found on the Adept Medical Website: www.adeptmedical.com/contact/serious-incident-reporting-form

The instructions provided above have been validated by the manufacturer of the medical device as being capable of preparing a medical device for reuse. It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process.

Visit the Adept Medical Website for detailed visual instructions related to this device. www.adeptmedical.com

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