# STARSYSTEM

AM0160 • AM0140 • AM0100

## INSTRUCTIONS FOR USE

Issue Date: 03/2026

Revision: H

**EC Representative**  
Qualrep Services B.V.  
Utrechtseweg 310 Bldg B42  
6812 AR Arnhem  
Netherlands

**Manufacturer**  
Adept Medical Ltd  
2-6 McDonald St,  
Morningside  
Auckland 1025, New Zealand

[adeptmedical.com](https://adeptmedical.com) |

ADEPT

Medical  
Supporting you

**Product Registration**

Register today to ensure your product is covered by warranty, and for easy access to product services and general support.

[www.adeptmedical.com/contact/product-registration](http://www.adeptmedical.com/contact/product-registration)

**Warranty**

**Service**


**Support**

**Contents**

Unboxing	<b>4</b>	
STARBoard	<b>4</b>	
Assembly	<b>4</b>	
Setup	<b>5</b>	
STARTable	<b>6</b>	
Setup	<b>6</b>	
STARSupport	<b>7</b>	
Assembly	<b>7</b>	
Setup	<b>7</b>	
Spare Parts	<b>11</b>	
Image Table Verification	<b>12</b>	
Essential Information	<b>14</b>	
ES   Spanish	<b>16</b>	NO   Norwegian <b>38</b>
FR   French	<b>18</b>	PL   Polish <b>40</b>
IT   Italian	<b>20</b>	PT   Portuguese <b>42</b>
DE   German	<b>22</b>	RO   Romanian <b>44</b>
NL   Dutch	<b>24</b>	RU   Russian <b>46</b>
HR   Croatian	<b>26</b>	SR   Serbian <b>48</b>
CS   Czech	<b>28</b>	SK   Slovak <b>50</b>
DA   Danish	<b>30</b>	SL   Slovenian <b>52</b>
FI   Finnish	<b>32</b>	SV   Swedish <b>54</b>
EL   Greek	<b>34</b>	TR   Turkish <b>56</b>
HU   Hungarian	<b>36</b>	

## Unboxing

	<b>A</b>	STARBoard x 1
	<b>B</b>	STARTable Shield x 1
	<b>C</b>	STARTable Shelf x 1
	<b>D</b>	STARSupport x 1
	<b>E</b>	Arm Support Pad x 1

AM0160 (STARSystem)

- A B C D E**

AM0140 (STARTable)

- B C**

AM0100 (STARBoard + STARSupport)



- A D E**

## STARBoard


## Assembly



Check friction joints and adjust if necessary


## Armrest Removal


## Setup


## Draping (Optional)

## Taping (Optional)




## Hand Positioning Options

STARTable

Setup




STARSupport



Assembly


Setup


Optional


5

Optional

7



Adjusting the STARSupport


1


6

2


Optional

STARSupport Detachment


Spare Parts


Product Codes

	STARSystem	AM0160
	STARBoard + STARSupport	AM0100
	STARTable	AM0140
	Arm Support Pad	AM0171
	STARSupport	AM0130

Component Codes


STARBoard


Maintenance instructions:  
[www.adeptmedical.com/imaging-table-accessories/starsystem/training?group=maintenance-instructions](http://www.adeptmedical.com/imaging-table-accessories/starsystem/training?group=maintenance-instructions)

**STARBoard Component Codes**


**STARTable**

**STARSupport**

**Image Table Verification**

The product has been risk assessed to operate within the below specifications

**Table Angulation:** +/- 15° lateral & longitudinal  
**Table Surface:** Flat

## Essential Information

### Information

#### Intended Purpose

STARBoard:

- To support a supine patient's arm, facilitating hyperextension for radial access.

STARSupport:

- To support the patient's left arm after radial artery access has been achieved.

STARTable:

- To provide an adjustable work surface.
- To provide scatter radiation protection for clinician.

#### Intended User/Training Requirement

Intended to be used by trained medical professionals.

#### Patient Target Group

Adult male or female patient requiring diagnostic imaging and image guided therapy.

#### Contraindications

- STARBoard and STARSupport: Not to be used with patients over 135 kg.
- STARTable: Shelf not to support weight in excess of 3 kg.
- Not for paediatric use.

#### Warnings and Cautions

- Ensure IFU is read prior to use.
- Ensure Preparatory Cleaning is conducted.
- Ensure Preparatory Product Check is conducted.
- Do not use with patients over 135 kg.
- Shelf does not support in excess of 3 kg.
- To provide restraint for restless patient (refer to page 5).
- Not for paediatric use.

#### Incident Reporting Guidelines

For product complaints and incidents, please complete form on [www.adeptmedical.com/contact/service-and-repairs](http://www.adeptmedical.com/contact/service-and-repairs)

#### Initial Checks & Cleaning

##### Cleaning

Refer to Disinfection Instructions.

##### Product Checks

- No exposed Carbon Fibre.
- No visible damage or sharp edges, e.g. cracks.
- STARBoard: Ensure all pivoting joints can rotate with the appropriate amount of friction and adjust if necessary (refer to page 4).
- STARBoard: Ensure correct placement (refer to page 5).
- STARBoard: Ensure the patient can access the Handle to grip (refer to page 5).
- STARTable: Ensure correct placement (refer to page 6).
- STARTable: Ensure Shelf is securely inserted (refer to page 6).
- STARSupport: Ensure that when locked, the device does not move.
- STARSupport: Ensure correct placement (refer to page 8).
- Arm Support Pad: Ensure no visible damage, or tearing.

#### Disinfection

##### Warnings

- Insufficient cleaning may compromise disinfection process and lead to transmission of pathogens.
- Do not use any disinfectant products or abrasive/corrosive agents which are not on the approved cleaners list.
- Always read manufacturer's instructions and consult the manufacturer's MSDS for cleaning and disinfectant products.
- Do not use any cleaning procedure which exceeds 65°C/149°F.
- Some disinfectants may cause slight discolouration to the soft blue material used on some components within the product range. This will not affect the strength and the product will remain fit for purpose.

##### Limitations on Processing

Discontinue use if:

- Any cracks, breakages or tears are present.
- Any exposed carbon fibre is visible.

#### Disinfection Instructions

##### Initial Treatment at the Point of Use

It is important to clean the product once it is removed from the packaging and after each procedure. Ensure all areas including joints, clips, sockets, brackets and levers are thoroughly cleaned after use to remove all contaminant build-up that may be present following a procedure.

##### Preparation Before Cleaning

- STARBoard: Remove STARBoard from beneath the mattress.
- Arm Support Pad: Remove Arm Support Pad from the STARBoard.
- STARSupport: Detach STARSupport from STARBoard and remove the swivel clamp.
- STARTable: Disassemble the Shelf and Shield.

[www.adeptmedical.com/imaging-table-accessories/starsystem/training?group=cleaning-instructions](http://www.adeptmedical.com/imaging-table-accessories/starsystem/training?group=cleaning-instructions)

##### Cleaning: Manual

- If any visible residue is present, rinse under water, but avoid submerging the device.
- A soft bristle brush can be used.
- Ensure the device is completely dried prior to disinfection.

#### Disinfection

Refer to the Approved Cleaners List for selection of an appropriate disinfecting agent.

STARBoard:

- With a cloth dampened by an approved disinfectant, wipe all surfaces, Armrest, Handle, Wrist Pivot, and adjust device to wipe entire Slider Arm.

STARSupport:

- With a cloth dampened by an approved disinfectant, wipe the cup and support. Then also wipe the lever, clamp, and Armrest.

STARTable:

- With a cloth dampened by an approved disinfectant, wipe all surfaces, brackets and spigots.

Arm Support Pad:

- With a cloth dampened by an approved disinfectant, wipe all surfaces.

#### Approved Cleaners List

##### Approved Disinfectant by Brands

###### Liquids:

ORION Laboratories - 70% Isopropyl Alcohol  
Jaychem Industries - 2% Chlorhexidine / 70% Alcohol  
Betadine - Povidone Iodine (7.5%)  
BODE Chemie GmbH - Kohrsolin FF  
STERIS - Coverage Spray TB  
CaviCide - Metrex Research

###### Dissolvable:

Du Pont - Rely+On Virkon Tablets

###### Surface Wipes:

PDI Sani-Cloth Bleach  
PDI Sani-Cloth AF3  
Reynard Health Supplies - Surface Disinfectant Wipes  
Metrex Research - CaviWipes  
Clinell Wipes - Universal (Green)  
Clorox Hydrogen Peroxide Wipes

##### Approved Active Disinfectant Ingredients

###### Acid Based:

≤10% Malic Acid CAS 6915-15-7  
≤6% Sulfamic acid CAS 5329-14-6

###### Alcohol Based:

≤5% 2-Butoxyethanol CAS 111-76-2  
≤10% Butyldiglycol CAS 112-34-5  
≤70% Isopropyl Alcohol (Propan-2-ol) CAS 67-63-0  
≤10% Tridecanol CAS 69011-36-5  
≤10% Alcohols, C12-14, ethoxylated CAS 68439-50-9  
≤70% Denatured Ethanol CAS 64-17-5

###### Sulphate Based:

≤0.1% PHMB CAS 27083-27-8  
≤55% Potassium Peroxymonosulfate CAS 70693-62-8  
≤3% Potassium Persulfate CAS 7727-21-1

###### Ammonium and Chloride Based:

≤5% Benzalkonium Chloride CAS 68424-85-1  
≤0.28% Benzethonium Chloride CAS 121-54-0  
≤10% Benzyl-C23-18-Alkyl-dimethyl Ammonium Chloride CAS 8001-54-5  
≤2% Chlorhexidine CAS 55-56-1  
≤10% Didecyl Dimethyl Ammonium chloride CAS 7173-51-5  
≤0.5% Quaternary Ammonium Compounds CAS 68956-79-6  
≤5% Benzyl-C 12-18 alkyl dimethyl ammonium chlorides CAS 63891-01-5

###### Others:

≤10% Glutaral CAS 111-30-8  
≤0.63% Sodium Hypochlorite CAS 7681-52-9  
≤7.5% Povidone Iodine CAS 25655-41-8  
0.5% Disodium Cocoampho Dipropionate CAS68604-71-7  
≤1.4% Hydrogen Peroxide CAS 7722-84-1

Access the latest Approved Cleaners list via [www.adeptmedical.com/resources](http://www.adeptmedical.com/resources)

#### Drying

All components should be dried thoroughly before use.

#### Maintenance, Inspection and Testing

All components:

- Ensure no exposed Carbon Fibre is visible.
- Ensure no visible damage or sharp edges, e.g. cracks.

STARBoard:

- Ensure all pivoting joints can rotate with the appropriate amount of friction and adjust if necessary.

STARSupport:

- Ensure when locked the Left Arm Support does not move.

#### Packaging

It is not required to package the device following disinfection.

#### Sterilisation

This device should not be subjected to sterilisation processes.

#### Storage

Once disinfection is complete and all components are dry the device should be stored in a dry environment.

#### Disposal

Used product is a biohazard, decontaminate according to instruction provided in the Disinfection Instructions and reach out to the manufacturer on [adeptmedical@adept.co.nz](mailto:adeptmedical@adept.co.nz) for material specification, if required. Dispose according to internal clinic policy taking into consideration local regulations.

#### Serious Incidents

Any serious incident which occurs in relation to the device should be reported immediately to:

- The Competent Authority of the applicable Member State.
- The Manufacturer using the Serious Incident Reporting Form found on the Adept Medical Website: [www.adeptmedical.com/contact/serious-incident-reporting-form](http://www.adeptmedical.com/contact/serious-incident-reporting-form)

The instructions provided above have been validated by the manufacturer of the medical device as being capable of preparing a medical device for reuse. It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process.

Visit the Adept Medical Website for detailed visual instructions related to this device. [www.adeptmedical.com](http://www.adeptmedical.com)

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[www.adeptmedical.com/imaging-table-accessories/starsystem/training?group=cleaning-instruetria](http://www.adeptmedical.com/imaging-table-accessories/starsystem/training?group=cleaning-instruetria)

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#### Liquidi:

ORION Laboratories - Alcol isopropilico 70%
Jaychem Industries - Cloroexidina 2%/Alcol 70%
Betadine - Iodopovidone (7,5%)
BODE Chemie GmbH - Kohrsolin FF
STERIS - Coverage Spray TB
CaviCide - Metrex Research

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