




CLINICIAN ARM SUPPORT

For supporting the clinician's forearm during procedures utilising endoscopic technology.

CLINICIAN ARM SUPPORT

For supporting the clinician's forearm during procedures utilising endoscopic technology.



The Adept Medical Clinician Arm Support offers an optimal solution for supporting the clinician's arm during surgeries where endoscopic technology is utilised.



Stable and Adjustable


The Clinician Arm Support is stable, height and angular adjustable, and can be tuned to suit individual preferences. The pivoting arms allow the forearm support to be locked rigid or left free to suit clinician and procedural requirements.



The support allows the clinician to manipulate the endoscope freely with low friction movement, providing a level of manoeuvrability. The Arm Rest comfortably supports the clinician's arm, offering an ergonomic solution to reduce strain and fatigue over the duration of a procedure.




With a quick and easy set-up, the Base Board placement allows for the support to be placed at any convenient position on the operating table, offering suitability for a range of endoscopic procedures. The cantilevered design frees up floorspace in the working environment, allowing clinicians to move easily.



Key Features and Benefits

- | Clinically Designed
- | Fits all operating tables
- | Quick set-up
- | Easy positioning
- | Radiolucent
- | Crafted from Carbon Fibre
- | Height and angular adjustable
- | Ergonomic

Specifications



 Clinician Arm Support
1.8 kg

Contact us for further details
adeptmedical@adept.co.nz

Applications

- Sinus Surgery (FESS)



Fields you can use this device



Ear, Nose
& Throat



Manufacturer

Adept Medical Ltd
2-6 McDonald St,
Morningside
Auckland 1025, New Zealand