



MR Safe


OVERHEAD ARM SUPPORT

Product Comparison

The Adept Medical Overhead Arm Support is now available in two sizes. It is the ideal imaging positioning solution to secure your patient in the arms overhead position to improve access and imaging to the patient's abdomen, removing unwanted image artefacts caused when the arms are left by their side, during C-arm, CT and MR.

Capable of holding one or both arms in the supine position, the Overhead Arm Support is a suitable arm positioning device for patients under general anaesthesia, or awake.

Issue Date: 02/2026



MR




NEW Overhead Arm Support MR Safe – Small
AM4100



MR

Overhead Arm Support MR Safe
AM4000

Product Specifications

Device	Overhead Arm Support MR Safe		Overhead Arm Support MR Safe – Small
Product Codes	AM4000		AM4100
Patient Weight	47 - 135 kg (103.5 - 297.5 lbs)		13 - 47 kg (28.5 - 103.5 lbs)
Dimensions	 *Height adjustable from 265 - 376 mm (10.5" - 15")		 **Height adjustable from 199 - 310 mm (8" - 12")
Weight	1.9 kg (4 lbs)		1.5 kg (3.5 lbs)
CT/MRI Bore Size	≥ 700 mm (27.5")		≥ 600 mm (23.5")
Imaging Table Compatibility	Table Width	450 - 700 mm (17.5" - 27.5")	
	Table Length	≥ 2000 mm (78.5")	
	Table Angulation	+/-15° lateral & longitudinal	
	Table Surface	Flat/Curved	
Aluminium Equivalence	 0.63 mm		

